



# Life Science




Cisa Group produces washing and sterilization equipment for the **healthcare** and **life science** industries.

Cisa **Life Science** equipment are designed for scientific laboratories, research centers, and pharmaceutical companies, active in the development and manufacturing of life-enhancing products.



**Cisa Life Science**  
info@cisagroup.it



Washing

Sterilization

Waste management

## Bridging Science and Life

Cisa **Life Science** supports scientists, researchers and engineers involved in the study and treatment of diseases.



Disease and drug discovery

Research



Clinical trials



Production



Therapy and pharma



# Applications



## Laboratories

- Washers  
Sterilizers
- BSL1/2
  - BSL3/4



## Research centers

- Sterilizers
- BSL1/2
  - BSL3/4
- Waste management



## Pharmaceutical


- Sterilizers




**Cisa Life Science**  
info@cisagroup.it

## Life Science products & solutions


Application	Laboratory	Research center	Pharma	External dimension 1P -2P	Chamber dimension		Biohazard configuration	
				Washer: measures are expressed in mm (W x H x D) Sterilizers: measures are expressed in mm (W x H)	Washer: measures are expressed in mm (W x H x D) and volume in L Sterilizers: measures are expressed in mm (W x H) and volume in L			
Washers	Small size mod. 104	●		780X2000x800-830	550X660x620	225		
	Medium size mod. 155	●		1200X2000x1036-1066	630X680x840	400		
	Large size mod. 305	●		1200X2000x1910	630X680x1680	800		
Sterilizers	P-3000LS	●		700x1850	322x322	71   101	✓	
	P-3600LS	●	●	903x1850	333x666	157   218	✓	
	P-420LS	●	●	903x1850	452x452	144   199   255	✓	
	P-450LS	●	●	1100x1850	452x660	207   296   379	✓	
	P-640LS	●	●	●	1424x1850	660x660	313   434   556   695   868	✓
	P-1000LS	●	●	●	2000x2100	660x1120	739   956   1178   1473   1843	✓
	P-1400LS		●	●	2000x2350	660x1490	982   1277   1572   1964   2455	✓
	P-1350LS		●	●	2800x2100	1050x1350	2268   2835   3544	✓
	P-2000LS		●	●	2800x2750	1050x2000	3358   4198   5247	✓
					Production capacity		Shredding	
Waste	WSD 100		●		Up to 100 Kg/H	1 shredder		
	WSD 200		●		Up to 200 Kg/H	1 shredder		
	WSD 400		●		Up to 400 Kg/H	2 shredders of 150/200 kg		
	WSD Sterily-Station WSM50		●		Up to 50 Kg/H	1 shredder		
	WSD Sterily-Station WSM70		●		Up to 70 Kg/H	1 shredder		




TRAINING



SERVICE



SOFTWARE



CONSUMABLES

## Related products